## INSTRUCTIONS

### HOW TO USE THE BUNSEN BURNER

**For KEMEX OUTFITS Nos. 2B and 3B**

This Outfit includes a Bunsen Burner and the necessary Rubber Tubing in place of a Spirit Lamp. This Burner can be used in all the experiments described in this Manual in which the Spirit Lamp is used as the source of heat.

The Bunsen Burner gives two types of flame. A yellow, luminous flame is obtained when the air hole at the base of the upright tube is closed by means of the rotating sleeve; a pale blue, non-luminous flame is produced when the air hole is open.

The non-luminous flame is hotter than the luminous one, and is used in nearly all

chemical experiments. This flame has an inner cone of unburned gas, and Test Tubes that are being heated should be held with their lower ends just above the tip of this cone. In using the Evaporating Stand the Ring should be fixed at the top of the Pillar, and the size of the flame adjusted by means of the gas tap so that the tip of the inner cone is just below the Gauze Square.

The luminous flame is used for experiments with the Blowpipe described on pages 2 and 3 of the No. 2-3 Manual. The gas tap should be partly closed in order to obtain a steady flame of medium size.

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